

# iFIX Head Mask Size Guide


Selecting the appropriate head mask size for your patient.

1




Measure the **distance between the lower chin and the top of the eye socket** with a medical tape measure.

2





Select the appropriate **head mask size** based on the closest dimension provided in the overview below.

3



Potential deviations will be **compensated** once the head mask is affixed to the selected anchor product.

## HEAD MASK SIZE GUIDE



Once you gain sufficient experience, you will select the right head mask easily without taking measurements. If you select the wrong size, dispose of the sample according to your local or hospital guidelines.